



## Lonigutamab Biosimilar - Anti-IGF1R mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1680
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1,Kappa
<b>SOURCE</b>	CAS 2362015-67-4

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Humanized
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Lonigutamab,IMMUNOGLOBULIN G1 (DE-C-TERMINAL-LYSINE), ANTI-(HUMAN INSULIN-LIKE GROWTH FACTOR I RECEPTOR) (HUMAN-MUS MUSCULUS MONOCLONAL HZ208F2-4 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL HZ208F2-4 .KAPPA.-CHAIN, DIMER, IMMUNOGLOBULIN

