



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Lorukafusp Biosimilar - Anti-Ganglioside GD2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1681
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Kappa
SOURCE	CAS 2131168-99-3

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Lorukafusp,ANTI-GANGLIOSIDE-GD2-ANTIBODY-INTERLEUKIN-2-FUSION-PROTEIN, APN301, EMD 273063, HU14.18-IL2, IMMUNOGLOBULIN G1-KAPPA, ANTI-(HUMAN GANGLIOSIDE GD2) CHIMERIC MONOCLONAL ANTIBODY, FUSED WITH INTERLEUKIN 2, GLYCOFORM ALFA, APN-301,EMD-273063,HU 14.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.