



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Melrilimab Biosimilar - Anti-IL1RL1 mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA1683 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG2,Kappa |
| SOURCE | CAS 2222865-46-3 |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Homo sapiens |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | PBS buffer PH7.5 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | Melrilimab,CNTO-7160,GSK3772847, GSK-3772847, IMMUNOGLOBULIN IGG2, ANTI-(HUMAN INTERLEUKIN 1 RECEPTOR-LIKE 1 PROTEIN) (HUMAN MONOCLONAL GSK3772847 .GAMMA.2-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL GSK3772847 .KAPPA.-CHAIN, DIMER,IL1RL1,anti-IL1RL1 |



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.