



Miromavimab Biosimilar - Anti-Rabies Virus Strain ERA GP Ectodomain Epitope G-II mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1687
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Kappa
SOURCE	CAS 2247163-73-9

PRODUCT INFORMATION

ORIGIN SPECIES	Mus musculus
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Miromavimab,IMMUNOGLOBULIN G1, ANTI-(RABIES VIRUS STRAIN EVELYN-ROCKITNIKI-ABELSETH GLYCOPROTEIN ECTODOMAIN EPI TOPE G-II) (MUS MUSCULUS MONOCLONAL M777-16-3 .GAMMA.1-CHAIN), DISULFIDE WITH MUS MUSCULUS MONOCLONAL M777-16-3 .KAPPA.-CHAIN, DIMER,Rabies Viru

