



Modakafusp Biosimilar - Anti-ADPRC mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1688
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4,Kappa
SOURCE	CAS 2254522-19-3

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Modakafusp,ANTI-CD38-TARGETED IGG4-ATTENUATED IFNA TAK-573, MODAKAFUSP ALFA, TEV-48573,ADPRC,anti-ADPRC

