



Obrindatamab Biosimilar - Anti-CD276;CD3E mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1692
CLONALITY	Monoclonal Antibody
ISOTYPE	Bispecific scFv with Domain Crossover,Kappa;Lambda
SOURCE	CAS 2069959-72-2

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Obrindatamab,DUAL-AFFINITY B7-H3/CD3-TARGETED PROTEIN MGD009, MGD-009,CD276;CD3E,anti-CD276;CD3E

