



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Opucolimab Biosimilar - Anti-PDL1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1697
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Lambda
SOURCE	CAS 2251771-79-4

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Opucolimab,IMMUNOGLOBULIN G1, ANTI-(HUMAN PROGRAMMED DEATH-LIGAND 1) (HUMAN MONOCLONAL HLX20 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL HLX20 .LAMBDA.-CHAIN, DIMER,PDL1,anti-PDL1



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.