



## Pavurutamab Biosimilar - Anti-TNFRSF17;CD3E mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1700
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	Bispecific scFv, Kappa; Lambda; Fc-Tag
<b>SOURCE</b>	CAS 2250292-39-6

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Humanized
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Pavurutamab,AMG701, AMG-701, ANTI-BCMA X ANTI-CD3 BITE (BISPECIFIC T CELL ENGAGER) ANTIBODY CONSTRUCT,TNFRSF17;CD3E,anti-TNFRSF17;CD3E

