



## Ponsegromab Biosimilar - Anti-GDF15 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1704
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1,Kappa
<b>SOURCE</b>	CAS 2368950-15-4

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Humanized
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Ponsegromab,IMMUNOGLOBULIN G1 (237-ALANINE,238-ALANINE,240-ALANINE), ANTI-(HUMAN GROWTH DIFFERENTIATION FACTOR 15) (HUMAN-MUS MUSCULUS MONOCLONAL PF-06946860 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL PF-06946860 .KAPPA.-CHAIN, DIMER,GD

