



Ragifilimab Biosimilar - Anti-TNFRSF18 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1708
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Kappa
SOURCE	CAS 2207590-51-8

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Ragifilimab,ANTI-GITR AGONISTIC MONOCLONAL ANTIBODY INCAGN01876,TNFRSF18,anti-TNFRSF18

