



Revdofilimab Biosimilar - Anti-TNFRSF4 mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA1710 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1,Kappa |
| SOURCE | CAS 2230138-89-1 |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Chimeric |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | PBS buffer PH7.5 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | Revdofilimab, ABBV368, HUMANIZED, RECOMBINANT MONOCLONAL ANTIBODY THAT BINDS TO CELL SURFACE EXPRESSED ,OX40 (TNF RECEPTOR SUPERFAMILY MEMBER 4), ANTI-OX40 AGONIST MONOCLONAL ANTIBODY ABBV-368, PR1628103, PR-1628103, TNFRSF4, anti-TNFRSF4 |

