



Sudubrilimab Biosimilar - Anti-CD274 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1717
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Kappa
SOURCE	CAS 2387417-06-1

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Sudubrilimab,IMMUNOGLOBULIN G1-KAPPA, ANTI-(HOMO SAPIENS CD274 (PROGRAMMED DEATH LIGAND 1, PDL1, PD-L1, B7 HOMOLOG 1, B7H1)), MONOCLONAL ANTIBODY,CD274,anti-CD274

