



Sugemalimab Biosimilar - Anti-PDL1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1718
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4,Lambda
SOURCE	CAS 2256084-03-2

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Sugemalimab,ANTI-PD-L1 MONOCLONAL ANTIBODY CS1001, WBP 3155, IMMUNOGLOBULIN G4 (228-PROLINE) (ANTI-(HUMAN PROGRAMMED DEATH- LIGAND 1)) (HUMAN MONOCLONAL CS1001 .GAMMA.4-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL CS1001 .LAMBDA.-CHAIN, DIMER,PDL1,anti-PDL1

