



Tarlatamab Biosimilar - Anti-DLL3;CD3E mAb - Research Grade

OVERVIEW

| | |
|--------------------------|------------------------------|
| PRODUCT REFERENCE | PX-TA1720 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | Bispecific scFv;Kappa;Lambda |
| SOURCE | CAS 2307488-83-9 |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Humanized |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | PBS buffer PH7.5 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | Tarlatamab,IMMUNOGLOBULIN SCFV-SCFV-SCFC, ANTI-(HOMO SAPIENS DLL3 (DELTA-LIKE LIGAND 3)) AND ANTI-(HOMO SAPIENS CD3E (CD3 EPSILON, LEU-4)), MONOCLONAL ANTIBODY SINGLE CHAIN (SCFV)2-SCFC, BISPECIFIC,DLL3;CD3E,anti-DLL3;CD3E |

