



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Tebentafusp Biosimilar - Anti-CD3E mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1721
CLONALITY	Monoclonal Antibody
ISOTYPE	Kappa, Strep + His-Tag
SOURCE	CAS 1874157-95-5

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Tebentafusp,IMCGP 100, T CELL RECEPTOR .ALPHA. CHAIN (SYNTHETIC HUMAN) FUSION PROTEIN WITH T CELL RECEPTOR .BETA. CHAIN (SYNTHETIC HUMAN) FUSION PROTEIN WITH IMMUNOGLOBULIN, ANTI-(HUMAN CD3 ANTIGEN) (SYNTHETIC SCFV FRAGMENT),CD3E,anti-CD3E



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.