



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Tirnovetmab Biosimilar - Anti-IL31 (Canine) mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1726
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2,Kappa
SOURCE	CAS 2364504-80-1

PRODUCT INFORMATION

ORIGIN SPECIES	Caninized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Tirnovetmab,IMMUNOGLOBULIN G-B, ANTI-(CANIS FAMILIARIS INTERLEUKIN 31) (CANIS FAMILIARIS-MUS MUSCULUS MONOCLONAL KIND-016 .GAMMA.B-CHAIN), DISULFIDE WITH CANIS FAMILIARIS-MUS MUSCULUS MONOCLONAL KIND-016 .KAPPA.-CHAIN, DIMER,IL31 (Canine),anti-IL31 (Canin



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.