



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Upifitamab Biosimilar - Anti-SLC34A2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1732
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Kappa
SOURCE	CAS 2254118-43-7

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Upifitamab,IMMUNOGLOBULIN G1 (DE-450-LYSINE), ANTI-(HUMAN SODIUM-DEPENDENT PHOSPHATE TRANSPORT PROTEIN 2B) (HUMAN-MUS MUSCULUS MONOCLONAL XMT-1535 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL XMT-1535 .KAPPA.-CHAIN, DIMER,,SLC34A2,anti-SL



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.