



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Vixtimotamab Biosimilar - Anti-CD3E;CD33 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1738
CLONALITY	Monoclonal Antibody
ISOTYPE	Bispecific Homodimer (VK-VH-VL'-VH', Tandem diabody);Kappa;Lambda
SOURCE	CAS 2243775-32-6

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Vixtimotamab,AMV564, MABFRAG HUMAN (FAB) ANTI P20138 (CD33_HUMAN)MABFRAG HUMAN (FAB) ANTI P07766 (CD3_HUMAN) (AMV564), CNTO 3953 TANDAB,T-652,CNTO-3953, AMV-564,CD3E;CD33,anti-CD3E;CD33



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.