



Vudalimab Biosimilar - Anti-CTLA4;PDCD1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1739
CLONALITY	Monoclonal Antibody
ISOTYPE	Bispecific ,IgG1,Kappa;Kappa
SOURCE	CAS 2329669-72-7

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Vudalimab,ANTI-PD-1/ANTI-CTLA-4 XMAB20717, IMMUNOGLOBULIN HALF-IG G1-KAPPA/SCFV-H-CH2-CH3, ANTI-(HOMO SAPIENS CTLA4 (CYTOTOXIC T-LYMPHOCYTE-ASSOCIATED PROTEIN 4, CD152)) AND ANTI-(HOMO SAPIENSPDCD1 (PROGRAMMED CELL DEATH 1, PD-1, PD1, CD279)), MONOCLONAL

