



Socazolimab Biosimilar - Anti-PDL1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1751
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-lambda2
SOURCE	CAS: 2305043-30-3

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens , Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Socazolimab,STI-1014, STI-A1014, ZKAB-001, ZKAB001,PDL1,anti-PDL1

