



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Simridarlimab Biosimilar - Anti-PDL1 & MER6 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1753
CLONALITY	Monoclonal Antibody
ISOTYPE	(G1_L-kappa)_VH-VH-h-CH2-CH3
SOURCE	CAS: 2378862-90-7

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens , Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Simridarlimab, IBI 322, IBI-322, IBI322, PDL1 & MER6, anti-PDL1 & MER6



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.