



Gresonitamab Biosimilar - Anti-CLDN18 & CD3e mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1758
CLONALITY	Monoclonal Antibody
ISOTYPE	(scFv-heavy-kappa)-(scFv-heavy-lambda)-scFc
SOURCE	CAS: 2413817-97-5

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Gresonitamab,AMG 910, AMG-910, AMG910,CLDN18 & CD3e,anti-CLDN18 & CD3e

