



Vepsitamab Biosimilar - Anti-MUC17 & CD3e mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1759
CLONALITY	Monoclonal Antibody
ISOTYPE	(scFv-heavy-lambda)-(scFv-heavy-lambda)-scFc
SOURCE	CAS: 2413453-53-7

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Vepsitamab,AMG 199, AMG-199, AMG199,MUC17 & CD3e,anti-MUC17 & CD3e

