



## Rocatinlimab Biosimilar - Anti-CD134 mAb - Research Grade

### OVERVIEW

---

|                          |                     |
|--------------------------|---------------------|
| <b>PRODUCT REFERENCE</b> | PX-TA1784           |
| <b>CLONALITY</b>         | Monoclonal Antibody |
| <b>ISOTYPE</b>           | IgG1-kappa          |
| <b>SOURCE</b>            | CAS: 2431972-52-8   |

### PRODUCT INFORMATION

---

|                           |   |
|---------------------------|---|
| <b>ORIGIN SPECIES</b>     | Homo sapiens  |
| <b>EXPRESSION SYSTEM</b>  | XtenCHO   |
| <b>BUFFER</b>             | PBS buffer PH7.5  |
| <b>STORAGE CONDITIONS</b> | Store at 4 °C for short term. Freeze and store at -80 °C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| <b>ALIAS/SYNONYMS</b>     | Rocatinlimab,KHK 4083, KHK-4083, KHK4083,CD134,anti-CD134   |

