



Manfidokimab Biosimilar - Anti-IL-4 receptor subunit alpha mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1788
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	CAS: 2428381-56-8

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Manfidokimab,AK 120, AK-120, AK120,IL-4 receptor subunit alpha,anti-IL-4 receptor subunit alpha

