



Upanovimab Biosimilar - Anti-Spike RBD mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1798
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS: 2500633-46-3

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Upanovimab,SCTA 01, SCTA-01, SCTA01,Spike RBD,anti-Spike RBD

