



Cibisatamab Biosimilar - Anti-CEACAM5&CD3E;CD3E mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1824
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1 Kappa; IgG1 Lambda
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Bispecific mAb with Domain Crossover
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Cibisatamab,,CEACAM5&CD3E;CD3E,anti-CEACAM5&CD3E;CD3E

