



MEDI3185 Biosimilar - Anti-Stromal cell-derived factor 1 receptor mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1898
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Stromal cell-derived factor 1 receptor,LCR1,HM89,LPS-associated protein 3,LAP-3,CXC-R4,CXCR-4,SDF-1 receptor,CD184,NPYRL,LESTR,Fusin,FB22,C-X-C chemokine receptor type 4,CXCR4,Lipopolysaccharide-associated protein 3,Leukocyte-derived seven transmembr

