



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## M802 Biosimilar - Anti-Bispecific mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1899
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>SOURCE</b>	

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Humanized
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	anti-Bispecific,p185erbB2,HER2,NGL,Tyrosine kinase-type cell surface receptor HER2,NEU,MLN 19,Proto-oncogene Neu,MLN19,ERBB2,Proto-oncogene c-ErbB- 2,CD340,Receptor tyrosine-protein kinase erbB-2,Metastatic lymph node gene 19 protein,T3E,T-cell surface ant



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**