



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Fabrazyme/Agalsidase Alfa Biosimilar - Anti-Alpha-D-galactosidase A mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1908
CLONALITY	Monoclonal Antibody
ISOTYPE	
SOURCE	CAS: 104138-64-9

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER PBS buffer PH7.5

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-Alpha-D-galactosidase A, Galactosylgalactosylglucosylceramidase
GLA, GLA, Melibiase, Agalsidase, Alpha-D-galactoside galactohydrolase, Alpha-galactosidase A,



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.