



VAK-694 Biosimilar - Anti-Binetrakin mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1911
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
SOURCE	VAK-694, NVP-VAK-694

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Binetrakin, B-cell stimulatory factor 1, Pitrakinra, Lymphocyte stimulatory factor 1, IL-4, BSF-1, IL4, Interleukin-4, VAK-694, NVP-VAK-694

