



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

CM101 Biosimilar - Anti-MPIF2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1912
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG
SOURCE	CM-101, CM101

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-MPIF2, Eosinophil chemotactic protein 2, Small-inducible cytokine A24, CCL24, Eotaxin-2, SCYA24, MPIF-2, Myeloid progenitor inhibitory factor 2, CK-beta-6, C-C motif chemokine 24, CM-101, CM101



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.