



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

SAR113244 Biosimilar - Anti-CXC-R5 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1915
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	SAR113244

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-CXC-R5,CD185,MDR-15,C-X-C chemokine receptor type 5,CXCR5,CXCR-5,MDR15,Monocyte-derived receptor 15,Burkitt lymphoma receptor 1,BLR1,SAR113244



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.