



BMS-936559/MDX-1105 Biosimilar - Anti-B7-H1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1916
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4
SOURCE	BMS-936559, MDX-1105

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-B7-H1, Programmed cell death 1 ligand 1, PDCD1 ligand 1, PDCD1L1, B7 homolog 1, PDCD1LG1, PDL1, hPD-L1, Programmed death ligand 1, B7H1, PD-L1, CD274, BMS-936559, MDX-1105

