



Rozipafusp alfa Biosimilar - Anti-B7RP1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1920
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2-kappa
SOURCE	CAS: 2143449-47-0

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-B7RP1,B7-related protein 1,CD275,ICOSL,ICOSLG,KIAA0653,B7 homolog 2,ICOS ligand,B7RP-1,B7H2,B7-H2,B7-like protein GI50,

