



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Zatuximab Biosimilar - Anti-HER1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1921
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	Futuximab, DS 992, Sym 004 (combination of modotuximab and futuximab), CAS: 1310460-86-6

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-HER1, Receptor tyrosine-protein kinase erbB-1, Proto-oncogene c-ErbB-1, Epidermal growth factor receptor, EGFR, ERBB1, ERBB, Futuximab, DS 992, Sym 004 (combination of modotuximab and futuximab)



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.