



Pogalizumab Biosimilar - Anti-ACT35 antigen mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1923
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	Vonlerolizumab, MOXR-0916, MOXR0916, RG-7888, RO7021608, CAS: 1638935-72-4

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-ACT35 antigen, TXGP1L, TAX transcriptionally-activated glycoprotein 1 receptor, Tumor necrosis factor receptor superfamily member 4, TNFRSF4, OX40L receptor, CD134, Vonlerolizumab, MOXR-0916, MOXR0916, RG-7888, RO7021608

