



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Toralizumab Biosimilar - Anti-sCD40L mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1925
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	E6040, IDEC-131, CAS: 252662-47-8

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-sCD40L,Tumor necrosis factor ligand superfamily member 5,T-cell antigen Gp39,CD40L,TNFSF5,CD154,CD40-L,TNF-related activation protein,CD40 ligand,CD40LG,TRAP,E6040, IDEC-131



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.