



## **$\alpha/\gamma$ Antibody (Kwar23) Biosimilar - Anti-BIT mAb - Research Grade**

### **OVERVIEW**

---

<b>PRODUCT REFERENCE</b>	PX-TA1930
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1 silence
<b>SOURCE</b>	

### **PRODUCT INFORMATION**

---

<b>ORIGIN SPECIES</b>	Homo sapiens
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	anti-BIT, Signal-regulatory protein alpha-1, Brain Ig-like molecule with tyrosine-based activation motifs, MFR, SIRP, Sirp-alpha-1, Sirp-alpha-2, PTPNS1, SIRPA, SHP substrate 1, CD172 antigen-like family member A, Sirp-alpha-3, Inhibitory receptor SHPS-1,

