



BHQ880 Biosimilar - Anti-DKK1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1933
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa(silence)
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Dickkopf-1, SK, hDkk-1, Dkk-1, Dickkopf-related protein 1, DKK1

