



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Vopikitug Biosimilar - Anti-IL2-RA mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1937
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS: 2733581-99-0

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-IL2-RA, IL-2R subunit alpha, TAC antigen, IL-2-RA, p55, IL2RA, IL-2 receptor subunit alpha, Interleukin-2 receptor subunit alpha, CD25



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.