



Patecibart Biosimilar - Anti-Endothelin-1 receptor mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1962
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	CAS: 2642374-04-5

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Endothelin-1 receptor, ETA, ET-A, ETRA, Endothelin receptor type A, EDNRA, hET-AR, ETA-R

