



Trabikibart Biosimilar - Anti-Cytokine receptor common subunit beta mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1974
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	CAS: 2643974-98-3

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Cytokine receptor common subunit beta, IL3RB, IL5RB, CDw131, CD131, CSF2RB, GM-CSF/IL-3/IL-5 receptor common beta subunit

