



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Obertamig Biosimilar - Anti-HLAG mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1980
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	half IG G1-lambda2/scFv-h-CH2-CH3
<b>SOURCE</b>	CAS: 2750126-32-8

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Human
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	0.01M PBS, pH 7.4
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	anti-HLAG, HLA-G, sHLA-G, HLA G antigen, MHC class I antigen G, HLA-6.0, HLA class I histocompatibility antigen, alpha chain G, T3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T-cell surface glycoprotein CD3 epsilon chain



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**