



Acimtamig Biosimilar - Anti-CD30L receptor mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1981
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG VH-V-kappa-VH'-V-lambda homodimer
SOURCE	CAS: 2738880-26-5

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-CD30L receptor, Ki-1 antigen, CD30, TNFRSF8, D1S166E, Tumor necrosis factor receptor superfamily member 8, Lymphocyte activation antigen CD30, CD16a, CD16A, FCGR3A, Fc-gamma RIII-alpha, FcRIIIa, IgG Fc receptor III-2, Low affinity immunoglobulin gamm

