



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Berahyaluronidase Alfa Biosimilar - Anti-Hyaluronidase PH-20 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1986
CLONALITY	Monoclonal Antibody
ISOTYPE	Human hyaluronidase isozyme-5/HYAL5 fragment engineered variant with the structural domain replaced by the corresponding region of human hyaluronidase isozyme-1/HYAL1.
SOURCE	CAS: 2636716-20-4

PRODUCT INFORMATION

ORIGIN SPECIES	
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Hyaluronidase PH-20, Sperm surface protein PH-20, HYAL3, Sperm adhesion molecule 1, PH20, SPAM1, Hyal-PH20, Hyaluronoglucosaminidase PH-20



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.