

Elritercept Biosimilar - Anti-ACVR2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1991
CLONALITY	Monoclonal Antibody
ISOTYPE	Fusion-[human activin receptor type 2-A (activin receptor type IIA, ACTR-IIA) extracellular domain (ECD) fragment (1-109) engineered varian]2-IgG1-Fc (Fragment constant)
SOURCE	CAS: 2727114-26-1

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-ACVR2, ACTRIIA, ACTR-IIA, Activin receptor type IIA, ACVR2A, Activin receptor type-2A

