



Acazicolcept Biosimilar - Anti-CD28 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA1999

CLONALITY

ISOTYPE Fusion - [ICOSLG N-terminal fragment (1-122)]₂ - IGHG1 Fc (Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-CD28, TP44, T-cell-specific surface glycoprotein CD28, CD278, ICOS, Inducible T-cell costimulator, Activation-inducible lymphocyte immunomediatory molecule, ALLIM

