



Alefacept Biosimilar - Anti-CD2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2000

CLONALITY

ISOTYPE Fusion - [CD58 (LFA3, lymphocyte function associated antigen 3)]₂ - IGHG1 Fc
(Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-Erythrocyte receptor, SRBC, T-cell surface antigen CD2, Rosette receptor, CD2,
LFA-3 receptor, T-cell surface antigen T11/Leu-5, LFA-2

