

Atacicept Biosimilar - Anti-BAFF and APRIL mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2002

CLONALITY

ISOTYPE Fusion - [TNFRSF13B (tumor necrosis factor receptor (TNFR) superfamily member 13B, TACI, transmembrane activator and CAML Interactor)]2 - IGHG1 Fc (Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-TNF- and APOL-related leukocyte expressed ligand 2, APRIL, CD256, Tumor necrosis factor ligand superfamily member 13, TNF-related death ligand 1, TALL-2, TNFSF13, ZTNF2, A proliferation-inducing ligand, TALL2, TRDL-1, BAFF, ZTNF4, BLYS, Tumor necrosi

